



Application

Purification of plasmids using FastGene® Plasmid Mini Kit

Product

FastGene® Plasmid Mini Kit (FG-90402)

Manufacturer

NIPPON Genetics EUROPE

The following data was kindly provided by Ken Honda, Department of Pharmacology, Yamaguchi University, Graduate School of Medicine, Japan. We tried transfection (luciferase activity) with a plasmid extracted with FastGene® Plasmid Mini Kit.

Experimental conditions

• Sample

E. coli strain: TOP10

Vector: pGL4 Vector 4.3kb (Insert size 0.5~2kb) * Insert: Promoter region

Input amount: LB medium 2.5 mL

• Procedure

1. Plasmid DNA-transformed TOP10 was cultured in LB medium
2. *E. coli* was recovered from the culture (2.5 mL) and plasmid DNA was purified with FastGene® Plasmid Mini Kit (elution buffer 50 µL)
3. DNA yield was measured with NanoDrop (Result ①)
4. Plasmids were digested by restriction enzymes. Presumed fragment pattern was confirmed by electrophoresis (Result ②)
5. After transfection, luciferase activity was measured (Result ③)
Detection device: Perkin Elmer 2030 ARVO X series multi label reader

Result


① DNA yield and purity

#	Elution buffer volume (µL)	yield (µg)	Nucleic Acid Conc. Unit	A260	A280	A260/A280	A260/A230
pGL4-prom-498	50.0	16.3072	332.8 ng/µL	6.656	3.511	1.90	2.20
pGL4-prom-998	50.0	14.2541	290.9 ng/µL	5.818	3.088	1.88	2.17
pGL4-prom-1489	50.0	16.2092	330.8 ng/µL	6.617	3.655	1.81	1.97
pGL4-prom-1994	50.0	23.4465	478.5 ng/µL	9.571	5.063	1.89	2.24

② Electrophoresis

The purified plasmid was used for restriction enzyme check


Conditions: 1 lane 10 µL (250 ng) used, TAE buffer, 0.7% TAE agarose, 100V · 30 min



The numerical value is the assumed base length (bp)
* Bands below 500 bp are not detected in gel image.

③ Measurement of luciferase assay

A plasmid with the confirmed fragment was transfected into mouse myoblast C2C12 cells and a reporter gene assay was performed: 0.2 µg of DNA was transfected into 10⁵ cells by Lipofectamine 2000. Stimulation of cells after 24 h, the cell lysate was recovered and the luciferase assay was measured.



A A promoter-inserted vector

B Vector only (Without insert)

After stimulation, luciferase assay activity was measured in all cells.



Customer comment

I think that purification degree and recovery rate is reproducible. The cost efficiency is also great. Personally, I think it will be easier to use, if you can pull out the column from the collection tube a bit more smoothly, but handling is not complicated overall, and it was possible to handle multile samples in a short time without stress.