



# PRODUCT INFORMATION SHEET

## IMMOBILIZED ENDOPROTEINASE GLU-C (V8) G3M

# P3102

Endoproteinase Glu-C (V8) from *Staphylococcus aureus*.

Endoproteinase Glu-C hydrolyzes peptide and ester linkages specifically at the carboxyl end of glutamic acid (-Glu/-X; in ammonium carbonate pH 7.8, or ammonium acetate pH 4.0, buffer A) or of glutamic and aspartic acid (-Glu/-X and -Asp/-X; in phosphate buffer pH 7.8, buffer B).

**G3m:** 25 µg (22 units) endoproteinase Glu-C per CR-column immobilized on dextran.

This CR-column cuts at least 12 µg tubulin or 5 µg BSA per application in phosphate buffer.

<b>Nr.7 Storage buffer:</b>	50 mM Tris/HCl at pH 7.5, 5 mM EDTA .
<b>Nr.31 Reaction buffer:</b>	25 mM NH <sub>4</sub> -acetate, pH 4.0 (see above)
<b>Nr.32 Washing buffer:</b>	25 mM NH <sub>4</sub> -acetate, pH 4.0, 1 M NaCl
<b>Nr.62 Reaction buffer:</b>	50 mM phosphate buffer, pH 7.8 (see above)
<b>Nr.63 Washing buffer:</b>	50 mM phosphate buffer, pH 7.8, 1 M NaCl

### Protocol

For more details see MoBiTec-CRC-Handbook.

#### 1. Dilute delivered buffers (at least 2 ml each) with sterile doubly distilled water.

For 1 application you need

1 ml 10x reaction buffer and 9 ml doubly distilled water

2 ml 5x washing buffer and 8 ml doubly distilled water

1 ml 10x storage buffer and 9 ml doubly distilled water

The substrate should be in reaction buffer

#### 2. Equilibrate the CR-column with 10 ml reaction buffer.

Fill 10 ml reaction buffer into a syringe, let the reaction buffer run through the column by gravity to the upper filter. In case the buffer runs very slowly, apply pressure by a syringe.

#### 3. Load substrate solution in reaction buffer.

Small volumes (< 70 µl): spin the CR-column 5 seconds in a benchtop centrifuge (2000 rpm are sufficient). Let the substrate solution enter the matrix material.

Larger volumes: Let the substrate solution run through the column.

Flow-rate: up to 70 µl/minute

Keep the substrate in the column for about 1 minute at room temperature. Higher turn-over is obtained when the substrate is applied to the column again or incubated for longer times.

#### 4. Elute the product solution.

Small volumes (<70 µl): centrifuge the product out of the column.

Larger volumes: Let the substrate run through the column and spin the residual solution out of the matrix

Notice: Molecules < 700 Dalton have to be eluted with 7 ml reaction buffer.

It does not harm the columns if they run dry.

#### 5. Wash the column with 10 ml washing buffer.

#### 6. Equilibrate the column with 10 ml storage buffer.

Store the column at 4°C. Never freeze a CR-column!