




Cleaver SCIENTIFIC

identiTAPES



55 m Core Adapters

IdentiTape Dispensers are designed to accommodate tapes of both core diameters. The central spindle takes the smaller 12.7 m rolls directly while by using the supplied adapters, the larger 55 m rolls, shown in the adjacent photograph, allows both sizes to be used side by side.


Cleaver IdentiTapes are oil-proof, waterproof and acid resistant tapes which will withstand temperature extremes of -20°C to over +121°C. They will adhere to surfaces such as glass, plastic, metal, paper, rubber and any other dry surface yet can be easily removed with no adhesive residue being left on the original surface. IdentiTapes can be written on with ball point pen, pencil, felt markers.

Self-adhesive IdentiTapes are designed to provide the perfect solution to a wide variety of labelling requirements

- IDENTIFY SAMPLES EASILY
- OIL, WATER AND ACID RESISTANT
- PAPER BASED
- WITHSTANDS EXTREMES IN TEMPERATURES, -20°C TO +121.1°C
- IDEAL FOR USE IN A FREEZER, REFRIGERATOR OR INCUBATOR
- NO ADHESIVE RESIDUE LEFT WHEN REMOVED
- SURFACE CAN BE WRITTEN ON WITH PEN, PENCIL OR LAB MARKER

WWW.CLEAVERSCIENTIFIC.COM

identiTAPES



Length	Width	WHITE	GREEN	YELLOW	BLUE	RED	ORANGE	PINK	VIOLET
12.7m	12.7mm	ST-12/W	ST-12/G	ST-12/Y	ST-12/B	ST-12/R	ST-12/O	ST-12/P	ST-12/V
12.7m	19.0mm	ST-34/W	ST-34/G	ST-34/Y	ST-34/B	ST-34/R	ST-34/O	ST-34/P	ST-34/V
12.7m	25.4mm	ST-10/W	ST-10/G	ST-10/Y	ST-10/B	ST-10/R	ST-10/O	ST-10/P	ST-10/V
55m	12.7mm	ST-120/W	ST-120/G	ST-120/Y	ST-120/B	ST-120/R	ST-120/O	ST-120/P	ST-120/V
55m	19.0mm	ST-340/W	ST-340/G	ST-340/Y	ST-340/B	ST-340/R	ST-340/O	ST-340/P	ST-340/V
55m	25.4mm	ST-100/W	ST-100/G	ST-100/Y	ST-100/B	ST-100/R	ST-100/O	ST-100/P	ST-100/V
TD1		Acrylic Tape Dispenser				CA55	Additional core adapters, pack/5		

CLEAVER SCIENTIFIC LTD

Unit 41, Somers Road Industrial Estate,
Rugby, CV22 7DH
United Kingdom

T_ +44 (0)1788 565300
E_ INFO@CLEAVERSCIENTIFIC.COM
W_ WWW.CLEAVERSCIENTIFIC.COM